

Product Datasheet

Human dLK1 ELISA Kit (orb439146)

Description	Human Delta Like 1 Homolog (dLK1) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P80370
Sensitivity	0.216ng/mL
Entrez	8788
Expiration Date	Please enquire.