


Product Datasheet

CD10 antibody (orb43915)

Description	Mouse Monoclonal to CD10.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	mouse NALM-6 leukemia pre-B cell line (tissue/cell preparation)
Target	CD10
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Uniprot ID	P08473
Entrez	4311
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry (surface staining) of hum...