

Product Datasheet

Human DNAJC12 ELISA Kit (orb439159)

Description	Human DnaJ/HSP40 Homolog Subfamily C, Member 12 (DNAJC12) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9UKB3
Sensitivity	0.057ng/mL
Entrez	56521
Expiration Date	Please enquire.