

Product Datasheet

Bovine DNASE1 ELISA Kit (orb439160)

Description	Bovine Deoxyribonuclease I (DNASE1) Sandwich ELISA Kit
Reactivity	Bovine
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P00639
Sensitivity	0.32ng/mL
Entrez	282217
Expiration Date	Please enquire.