

Product Datasheet

Human DNMT3B ELISA Kit (orb439171)

Description	Human DNA Methyltransferase 3B (DNMT3B) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9UBC3
Sensitivity	0.054ng/mL
Entrez	1789
Expiration Date	Please enquire.