

Product Datasheet

Human DPEP1 ELISA Kit (orb439177)

Description	Human Dipeptidase 1, Renal (DPEP1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P16444
Sensitivity	0.050ng/mL
Entrez	1800
Expiration Date	Please enquire.