

Product Datasheet

Human DPP4 ELISA Kit (orb439181)

Description	Human Dipeptidyl Peptidase IV (DPP4) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P27487
Sensitivity	14.8pg/mL
Entrez	1803
Expiration Date	Please enquire.