

Product Datasheet

Mouse DPP9 ELISA Kit (orb439187)

Description	Mouse Dipeptidyl Peptidase 9 (DPP9) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8BVG4
Sensitivity	25pg/mL
Entrez	224897
Expiration Date	Please enquire.