

Product Datasheet

Human DRD1 ELISA Kit (orb439193)

Description	Human Dopamine Receptor D1 (DRD1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P21728
Sensitivity	0.109ng/mL
Entrez	1812
Expiration Date	Please enquire.