

Product Datasheet

Human DRD2 ELISA Kit (orb439195)

Description	Human Dopamine Receptor D2 (DRD2) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P14416
Sensitivity	0.129ng/mL
Entrez	1813
Expiration Date	Please enquire.