

Product Datasheet

Mouse DSC1 ELISA Kit (orb439200)

Description	Mouse Desmocollin 1 (DSC1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P55849
Sensitivity	0.051ng/mL
Entrez	13505
Expiration Date	Please enquire.