

Product Datasheet

Human DUSP1 ELISA Kit (orb439215)

Description	Human Dual Specificity Phosphatase 1 (DUSP1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P28562
Sensitivity	12.9pg/mL
Entrez	1843
Expiration Date	Please enquire.