

Product Datasheet

Mouse DUSP5 ELISA Kit (orb439218)

Description	Mouse Dual Specificity Phosphatase 5 (DUSP5) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q1HL35
Sensitivity	0.055ng/mL
Entrez	240672
Expiration Date	Please enquire.