

Product Datasheet

Human DYNC1H1 ELISA Kit (orb439221)

Description	Human Dynein, Cytoplasmic 1, Heavy Chain 1 (DYNC1H1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q14204
Sensitivity	0.113ng/mL
Entrez	1778
Expiration Date	Please enquire.