

Product Datasheet

Human DYRK1A ELISA Kit (orb439222)

Description	Human Dual Specificity Tyrosine Phosphorylation Regulated Kinase 1A (DYRK1A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q13627
Sensitivity	0.117ng/mL
Entrez	1859
Expiration Date	Please enquire.