

Product Datasheet

Monkey ECF ELISA Kit (orb439247)

Description	Monkey Eosinophil Chemotactic Factor (ECF) Sandwich ELISA Kit
Reactivity	Monkey
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8MIT7
Sensitivity	7.3pg/mL
Entrez	574218
Expiration Date	Please enquire.