

Product Datasheet

Human ECRG4 ELISA Kit (orb439250)

Description	Human Esophageal Cancer Related Gene 4 (ECRG4) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q99LS0
Sensitivity	30pg/mL
Entrez	78896
Expiration Date	Please enquire.