

Product Datasheet

Mouse EGLN1 ELISA Kit (orb439264)

Description	Mouse Egl Nine Homolog 1 (EGLN1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q91YE3
Sensitivity	0.108ng/mL
Entrez	112405
Expiration Date	Please enquire.