

Product Datasheet

Mouse EIF2aK3 ELISA Kit (orb439267)

Description	Mouse Eukaryotic Translation Initiation Factor 2 Alpha Kinase 3 (EIF2aK3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Z2B5
Sensitivity	0.108ng/mL
Entrez	13666
Expiration Date	Please enquire.