

Product Datasheet

Human EIF4A2 ELISA Kit (orb439268)

Description	Human Eukaryotic Translation Initiation Factor 4A2 (EIF4A2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q14240
Sensitivity	0.059ng/mL
Entrez	1974
Expiration Date	Please enquire.