

Product Datasheet

Human EIF6 ELISA Kit (orb439269)

Description	Human Eukaryotic Translation Initiation Factor 6 (EIF6) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P56537
Sensitivity	0.123ng/mL
Entrez	3692
Expiration Date	Please enquire.