

Product Datasheet

Human ELA1 ELISA Kit (orb439270)

Description	Human Elastase 1, Pancreatic (ELA1) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9UNI1
Sensitivity	23.6pg/mL
Entrez	1990
Expiration Date	Please enquire.