

Product Datasheet

Mouse ELA1 ELISA Kit (orb439271)

Description	Mouse Elastase 1, Pancreatic (ELA1) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q91X79
Sensitivity	26.5pg/mL
Entrez	109901
Expiration Date	Please enquire.