

Product Datasheet

Human ELN ELISA Kit (orb439285)

Description	Human Elastin (ELN) Sandwich ELISA Kit
Reactivity	Human
Range	0.468-30ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P15502
Sensitivity	0.21ng/mL
Entrez	2006
Expiration Date	Please enquire.