

Product Datasheet

Rat ELOVL1 ELISA Kit (orb439289)

Description	Rat Elongation Of Very Long Chain Fatty Acids Like Protein 1 (ELOVL1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q5U2Z8
Sensitivity	0.054ng/mL
Entrez	679532
Expiration Date	Please enquire.