

Product Datasheet

Human EMILIN1 ELISA Kit (orb439290)

Description	Human Elastin Microfibril Interface Located Protein 1 (EMILIN1) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9Y6C2
Sensitivity	0.64ng/mL
Entrez	11117
Expiration Date	Please enquire.