

Product Datasheet

Rat EML2 ELISA Kit (orb439294)

Description	Rat Echinoderm Microtubule Associated Protein Like Protein 2 (EML2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q6P6T4
Sensitivity	0.126ng/mL
Entrez	192360
Expiration Date	Please enquire.