

## Product Datasheet

### Human EMMPRIN ELISA Kit (orb439295)

|                          |                                                                                                |
|--------------------------|------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Extracellular Matrix<br>Metalloproteinase Inducer (EMMPRIN)<br>Sandwich ELISA Kit        |
| <b>Reactivity</b>        | Human                                                                                          |
| <b>Range</b>             | 31.25-2000pg/mL                                                                                |
| <b>Form/Appearance</b>   | Traditional                                                                                    |
| <b>Storage</b>           | All reagents should be kept according to<br>the instructions on the packaging upon<br>receipt. |
| <b>Note</b>              | For research use only                                                                          |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell<br>lysates and other biological fluids.                |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell<br>lysates and other biological fluids.                |
| <b>Uniprot ID</b>        | <b>P35613</b>                                                                                  |
| <b>Sensitivity</b>       | 11.4pg/mL                                                                                      |
| <b>Entrez</b>            | <b>682</b>                                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                                                |