

Product Datasheet

Human EMP3 ELISA Kit (orb439297)

Description	Human Epithelial Membrane Protein 3 (EMP3) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P54852
Sensitivity	0.099ng/mL
Entrez	2014
Expiration Date	Please enquire.