

Product Datasheet

Human EN2 ELISA Kit (orb439299)

Description	Human Engrailed Homeobox Protein 2 (EN2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P19622
Sensitivity	0.058ng/mL
Entrez	2020
Expiration Date	Please enquire.