

## Product Datasheet

### Mouse ENK ELISA Kit (orb439310)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Enkephalin (ENK) competitive inhibition ELISA Kit                                  |
| <b>Reactivity</b>        | Mouse                                                                                    |
| <b>Range</b>             | 93.75-6000pg/mL                                                                          |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma and other biological fluids.                                               |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.                                               |
| <b>Sensitivity</b>       | 25.35pg/mL                                                                               |
| <b>Expiration Date</b>   | Please enquire.                                                                          |