

Product Datasheet

Human ENPP2 ELISA Kit (orb439313)

Description	Human Ectonucleotide Pyrophosphatase/Phosphodiesterase 2 (ENPP2) Sandwich ELISA Kit
Reactivity	Human
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q13822
Sensitivity	32.28pg/mL
Entrez	5168
Expiration Date	Please enquire.