

Product Datasheet

Mouse ENSa ELISA Kit (orb439314)

Description	Mouse Endosulfine Alpha (ENSa) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P60840
Sensitivity	0.065ng/mL
Entrez	56205
Expiration Date	Please enquire.