

Product Datasheet

Human EP300 ELISA Kit (orb439317)

Description	Human E1A Binding Protein P300 (EP300) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q09472
Sensitivity	0.071ng/mL
Entrez	2033
Expiration Date	Please enquire.