

Product Datasheet

Rat EPO ELISA Kit (orb439335)

Description	Rat Erythropoietin (EPO) Sandwich ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P29676
Sensitivity	24.2pg/mL
Entrez	24335
Expiration Date	Please enquire.