

Product Datasheet

Human ERa ELISA Kit (orb439344)

Description	Human Estrogen Receptor Alpha (ERa) Sandwich ELISA Kit
Reactivity	Human
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P03372
Sensitivity	0.255ng/mL
Entrez	2099
Expiration Date	Please enquire.