

Product Datasheet

Mouse ERa ELISA Kit (orb439345)

Description	Mouse Estrogen Receptor Alpha (ERa) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P19785
Sensitivity	0.146ng/mL
Entrez	13982
Expiration Date	Please enquire.