

Product Datasheet

Rat ERa ELISA Kit (orb439347)

Description	Rat Estrogen Receptor Alpha (ERa) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P06211
Sensitivity	0.117ng/mL
Entrez	24890
Expiration Date	Please enquire.