

Product Datasheet

Human ErbB4 ELISA Kit (orb439350)

Description	Human V-Erb A Erythroblastic Leukemia Viral Oncogene Homolog 4 (ErbB4) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q15303
Sensitivity	0.062ng/mL
Entrez	2066
Expiration Date	Please enquire.