

Product Datasheet

Human EZH1 ELISA Kit (orb439382)

Description	Human Enhancer Of Zeste Homolog 1 (EZH1) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q92800
Sensitivity	0.113ng/mL
Entrez	2145
Expiration Date	Please enquire.