

Product Datasheet

Human F10 ELISA Kit (orb439389)

Description	Human Coagulation Factor X (F10) Sandwich ELISA Kit
Reactivity	Human
Range	93.75-6000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P00742
Sensitivity	39.5pg/mL
Entrez	2159
Expiration Date	Please enquire.