

Product Datasheet

Human F11 ELISA Kit (orb439392)

Description	Human Coagulation Factor XI (F11) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	P03951
Sensitivity	0.127ng/mL
Entrez	2160
Expiration Date	Please enquire.