

Product Datasheet

Mouse F13B ELISA Kit (orb439398)

Description	Mouse Coagulation Factor XIII B Polypeptide (F13B) Sandwich ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q07968
Sensitivity	1.02ng/mL
Entrez	14060
Expiration Date	Please enquire.