

Product Datasheet

Human Coagulation Factor II ELISA Kit (orb439402)

Description	Human Coagulation Factor II (F2) Sandwich ELISA Kit
Reactivity	Human
Range	9.375-600ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum.
Sample Types	serum.
Uniprot ID	P00734
Sensitivity	3.0ng/mL
Entrez	2147
Expiration Date	Please enquire.