

Product Datasheet

Human Coagulation Factor V ELISA Kit (orb439405)

Description	Human Coagulation Factor V (F5) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma and other biological fluids.
Uniprot ID	P12259
Sensitivity	23.9pg/mL
Entrez	2153
Expiration Date	Please enquire.