

Product Datasheet

Mouse Coagulation Factor V ELISA Kit (orb439406)

Description	Mouse Coagulation Factor V (F5) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma and other biological fluids.
Uniprot ID	088783
Sensitivity	0.145ng/mL
Entrez	14067
Expiration Date	Please enquire.