

Product Datasheet

Human Coagulation Factor VII ELISA Kit (orb439409)

Description	Human Coagulation Factor VII (F7) Sandwich ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma and other biological fluids.
Uniprot ID	P08709
Sensitivity	19.56pg/mL
Entrez	2155
Expiration Date	Please enquire.