

Product Datasheet

Mouse Coagulation Factor VII ELISA Kit (orb439410)

Description	Mouse Coagulation Factor VII (F7) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma and other biological fluids.
Uniprot ID	P70375
Sensitivity	37pg/mL
Entrez	14068
Expiration Date	Please enquire.