

Product Datasheet

Canine Coagulation Factor VIII ELISA Kit (orb439412)

Description	Canine Coagulation Factor VIII (F8) Sandwich ELISA Kit
Reactivity	Canine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma and other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	O18806
Sensitivity	1.15ng/mL
Entrez	403875
Expiration Date	Please enquire.