

Product Datasheet

Human Coagulation Factor VIII ELISA Kit (orb439413)

Description	Human Coagulation Factor VIII (F8) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma and other biological fluids.
Uniprot ID	P00451
Sensitivity	0.055ng/mL
Entrez	2157
Expiration Date	Please enquire.